

The logo consists of a grid of small white squares on the left, followed by the word "Microtrac" in a bold, italicized, white sans-serif font.

Microtrac

Total Solutions in Particle Characterization

Nanotrak Wave

**Particle Size, Zeta Potential
and Molecular Weight Analyzer**

Introducing **Nanotracs Wave**

Microtrac has been a pioneer in particle sizing technology for over 30 years! The Nanotracs Wave is the latest generation of sub-micron Particle Size and Zeta Potential analyzers from Microtrac! The enhanced design on the Nanotracs Wave features faster measurements, smaller particle size capability (to 0.8 nm), higher precision, higher accuracy and advanced software capabilities in a small, robust dynamic light scattering instrument with no moving components. Nanotracs Wave complies with the ISO 13321 standard for Dynamic Light Scattering.

The main features of the Nanotracs Wave are

Technology: Dynamic Light Scattering incorporating the patented Controlled Reference Method for advanced power spectrum analysis of Doppler shifts under Brownian Motion.

Range: Measurement capability from 0.8 to 6500 nanometers. Zeta Potential from -125 to +125mV
Sample Size: Typically less than 3ml in standard cell, or 0.2ml with small cell option.

Rapid Analysis: 15 to 30 second analysis times in most cases.

Analysis Simplicity: No "A Priori" or advance knowledge of the sample is required. Nanotracs will accurately report monomodal, multimodal, broad and narrow distributions without the need to select special calculation algorithms.

Accuracy: Measures using Mie scattering calculations for spherical particles and proprietary Mie calculations for non-spherical particles. **This feature is unique to Microtrac.**

Repeatability: Better than 1% from instrument to instrument for 100nm polystyrene.

Traceability: Particle Size measurements traceable to NIST standards.

Concentration: Capable to measure up to 40% by solids in some cases. This feature reduces the need to dilute samples which may affect particle size or colloid stability. Low concentration measurements can be as low as 0.1ppm for 200nm polystyrene.

Temperature Control: Precise temperature control using a Peltier control device.


Versatility: Nanotracs Wave probe can be external for use in Dip-N-Run or in-line applications. Nanotracs Wave is also compatible with the **Zetrator** titration device.

Security/Validation: Nanotracs Wave FLEX software is compatible with FDA 21 CFR Part 11 protocols. Full IQ/OQ/PQ validation documentation is available for Nanotracs Wave.


How Nanotracs Wave works

Particles suspended in a dispersing fluid are subject to random collisions with the thermally excited molecules of the dispersing fluid resulting in Brownian motion. The velocity and direction of the resulting motion are random but the velocity distribution of a large number of mono-sized particles averaged over a long period will approach a known functional form, in this case the size distribution of the particles.

In the Nanotracs Wave, light from a laser diode is coupled to the sample through an optical beam splitter in the Nanotracs probe assembly. The interface between the sample and the probe is a sapphire window at the probe tip. The sapphire window has two functions! Firstly, it reflects the original laser back through the beam splitter to a photodetector. This signal which has the same frequency as the original laser acts as a reference signal for detection, offering Heterodyne detection. Secondly, the laser passes through the sapphire window and is scattered by the particles which are in suspension but moving under Brownian motion. The laser is frequency shifted according to the Doppler effect relative to the velocity of the particle. Light is scattered in all directions including 180 degrees backwards. This scattered, frequency shifted light is transmitted through the sapphire window to the optical splitter in the probe to the photodetector. These signals of various frequencies combine with the reflected signal of un-shifted frequency (Controlled Reference) to generate a wide spectrum of heterodyne difference frequencies. The power spectrum of the interference signal is calculated with dedicated highspeed FFT (Fast Fourier Transform) digital signal processor hardware. The power spectrum is then inverted to give the particle size distribution.



The Nanotracs incorporates a highly accurate temperature sensor in the sample cell. By describing the fluid temperature and viscosity characteristics in the Nanotracs algorithm, these parameters can be included in determining accurate particle size distributions. Also, because the laser light needs only to penetrate approximately 100 microns into the sample to generate a power spectrum, the Nanotracs can accurately determine particle size distributions at significantly higher concentrations than other methods.



Zeta Potential Measurement

The Nanotractor Wave also has the capability to measure Zeta Potential using a specific sample cell for size as well as charge measurements, Zetapotential measurement capability is provided by additional hardware in the form of :

1. A voltage source with programmable amplitude and wave form.
2. An insulating sample cell, with optical probes opposed by electrodes.
3. Optical probes, with sample interface window consisting of typical sapphire, but with specialized metallic and semiconductor optical coatings applied.

Coated-window optical probes are paired with their opposite electrodes in an insulating sample cell. Excitation of the cell from the voltage source is applied between the optical probes and their electrodes, creating an electric field. Particle motion is analyzed while under the influence of the field. Particle size distribution is determined from the velocity distribution of particles suspended in a dispersing medium, using the principles of dynamic light scattering. The Nanotractor Wave analyzer measures the additional velocity imparted to the charged particles when placed in an electric field. Particle electrophoretic mobility is calculated from this additional velocity component. Zeta potential is calculated from mobility using accepted relationships between mobility and zeta potential. The relationship between zetapotential and mobility is given by the Smoluchowski equation:

$$\zeta = \mu\eta/\epsilon$$

where

ζ = zeta potential,

μ = mobility,

η = viscosity,

ϵ = dielectric constant

for water at 25degC,

$$\text{Zeta potential(mV)} = 12.8 \times \text{Mobility}(\mu/\text{sec}/\text{volt}/\text{cm})$$

Nanotractor Wave Sample Modules

The Nanotractor Wave provides the user with a versatile selection of Sample Modules depending on measurement requirements. Users can easily change from size measurement only to size and zeta potential measurement by changing to the relevant sample module. For users requiring temperature control there is also a module with a controlled Peltier device.


The sample is introduced to a removable teflon sample cell. The sample cell and probe tip can be sterilized for users in the biotechnology sector.

Flow through cells are available for interfacing with the Microtractor Zetrator facilitating comprehensive forward and backward Titration sequences for product stability applications.

Molecular Weight Measurement

The Nanotractor Wave measures Molecular Weight using the Hydrodynamic sample properties of the material and also through the generation of Debye Plots. From the Debye plot 2nd Virial coefficient can also be generated. The size of suspended molecules is determined by Dynamic Light Scattering technology, the foundation of the Nanotractor Wave DLS Analyzer measurement. The DLS technology measures the Brownian motion of the molecules and determines the molecular size through the application of the Brownian motion principles. The size is referred to as the hydrodynamic size.

The measurement of Molecular Weight ('MW') using the Debye plot depends upon an accurate measurement of the light that is scattered by a molecular suspension of known concentration, C. The ratio of total scattered-light to total incident-light is the Rayleigh ratio, R. The Debye MW expression relates the MW to the Rayleigh ratio and to the sample concentration:



For Microtractor DLS Analyzers, the angle of incident-to-scattered light is 180° (backscatter). The Rayleigh Ratio is calculated from the DLS Analyzer's measurement of the backscattered light intensity.



Nanotracs Wave: Technical Data

Particle Size	0.08 to 6500 nanometers (displayed in nanometers, microns or angstroms)	
Zeta Potential	-125 to +125 mV, Mobility -10 to +10 microseconds/volt/cm, size range 10nm to 20 microns	
Molecular Weight	1kDa to 20 MDa, Hydrodynamic, Debye Plot with 2nd Virial co-efficient calculation	
Sample Volume	Size only - 50 ul to 500 ul, with micro insert min 10 to 20 ul. Zeta 0.7 to 3 ml.	
Measurement Angle	180 Degrees	
Repeatability	1% or better for 100nm Polystyrene	
Concentration Limits	From ppb to 40% by volume in certain conditions.	
Laser	Laser Diode, 780 nm , 3mW Nominal, Class IIIB, no alignment required.	
Temperature Control	5 to 90 Degrees Centigrade +/- 0.2 Degrees using Peltier Device	
Sample Cell Modules	Nano, Zeta, Micro insert for Nano Cell, Flow thru cell with or without cap.	
Chemical Compatibility	Aqueous, polar and non polar solvents, organic solvents. Surfaces-Stainless steel, sapphire and Teflon	
pH Range	2 to 12 pH	
Conductivity Range	1 to 200 mS/cm	
Environmental	Operating Temperature	10 to 50 Degrees C
	Humidity	up to 90% non condensing
Dimensions	15L x 14W x 13H inches (38.1 x 35.5 x 33 cms)	
Weight	30 Lbs, approx 15 Kg	
Electrical	90 to 240 VAC, 47 to 63 HZ, 5 Amps	

Contact Details

For more information on Nanotracs Wave as well as other Microtracs products contact Microtracs Inc at (+1) 727 507 9770 or contact your local Microtracs Representative or log on to our website at www.microtracs.com