

CAP8[®]

Semi-Automatic
Production Scale
Capsule Filling
Machine



CAPSUGEL[®]
Quality
People and Products Working Together™

The cost-effective, low maintenance leader

CAP8® Capsule Filling Machine

Our easy to operate CAP8® is a low-maintenance and cost-effective filling machine that offers semi-automatic operation and reliable performance. The CAP8 is designed to allow you to scale up to production-level volumes without the need for substantial investment.

Specifications

Production Capacity*

Capsule Size	Speed per hour
5.....	17,000–25,000
4-3.....	20,000–29,000
2-1.....	17,000–25,000
0-0el.....	14,000–22,000
00-00el**.....	9,000–15,000
000**.....	7,000–11,000

AAA**, AA**, AAel**, A, B and C: Contact Capsugel's Technical Support Group for details.

* Rates depend on operator's skill level and density and type of material being filled.
** Conversion kit required.



Dimensions

136 cm x 77 cm x 191 cm high

Weight

440 kg (970 lbs)

Power Requirements

200-240 VAC 3-Phase on a 20-amp circuit;
WYO or Delta Power

Other

All orders include 1 set of change parts and choice of auger (Spade, One-Way, Two-Way, or Three-Way).

To bring your production on line immediately, we offer an installation package including two days of technical training and up to three days of operator training.



User-friendly Features

A wide array of innovative features ensures years of efficient filling and reliability.

- Table-top pneumatic closing station for easy integration of capsule finishing equipment
- Powder hopper is removable for cleaning
- Easy auger adjustment to change fill weights
- Easy closed capsule length adjustment
- Quick and accurate change part set-up
- “Toolless” changeover of rectifier parts
- Easy speed selection
- Large powder recovery station
- Robust electronic controls
- Easy to remove rotary table and ring carrier assemblies

GMP Features

The CAP8’s advanced technology and components ensure cleanliness and protect product integrity.

- All contact surfaces are made with FDA approved materials
- Stainless steel ring carrier, rotary table and powder hopper
- The internal components of the cabinet are protected against powder infiltration
- Easy to clean

Safety Features

The unique safety features of the CAP8 protect operators and mechanics.

- Remote E-stop for operator safety
- Guards and proximity switch on closing station
- Rectifier guard interlocks
- CE compliant



Optional Features

- Casters
- Fill ring rest tray
- Free-standing vibrator

Support

Capsugel’s highly experienced technical service team is always ready to assist customers. Offering hands-on support for shop-floor operators and managers, they can optimize production performance and help prevent problems before they arise. They will also work closely with your key production and maintenance people to identify areas of improvement and to develop machine-specific action plans.

For more information on the CAP8® Capsule Filling Machine and how it can help improve your efficiency, contact your Capsugel sales representative (see list of worldwide locations on back) or visit us on the web at www.capsugel.com.



AMERICAS

Greenwood, South Carolina
 Capsugel
 535 North Emerald Road
 Greenwood, SC 29646
 USA
 Tel: 888-783-6361 or
 864-223-2270
 Fax: 888-783-6360

Sales Office México City, México
 Pfizer Consumer Health Care
 México, S. de R.L. de C.V.
 Div. Capsugel
 Paseo de los Tamarindos #40
 Col. Bosques de las Lomas
 05120 México, D.F. México
 Tel: (5255) 9177 3120
 Fax: (5255) 5257 5174

Sales Office Rio de Janeiro, Brazil
 Laboratórios Pfizer Ltda.
 Div. Capsugel
 Avenida Das Américas #8445
 Salas 1202 to 1207
 Barra da Tijuca
 CEP: 22793-081
 Rio de Janeiro, RJ, Brazil
 Tel: +55-21-2432 5266
 Fax: +55-21-2432-5270

EUROPE

Bornem, Belgium
 Capsugel Bornem
 Rijksweg 11
 B-2880 Bornem
 Belgium
 Tel: 32-(0)3-890-05-11
 Fax: 32-(0)3-889-26-22

Colmar, France
 Capsugel Colmar
 10, Rue Timken
 F-68027 Colmar, Cedex
 France
 Tel: 33-(0)3-89-20-57-09
 Fax: 33-(0)3-89-41-48-11

Ploërmel, France
 Capsugel
 Z.I. Camagnon - BP 320
 56 803 Ploërmel, Cedex
 France
 Tel: 33-(0)2-97-72-10-00
 Fax: 33-(0)2-97-72-15-00

Cambridge, UK
 Capsugel
 Cambridge Research Centre
 Compass House, Vision Park
 Chivers Way, Histon, Cambridge CB4 9ZR
 United Kingdom
 Tel: 44-(0)-1223-598480
 Fax: 44-(0)-1223-598070

ASIA

Ayutthaya, Thailand
 Capsugel (Thailand) Co., Ltd.
 1/82 Moo 5, Rojana Industrial Park
 Rojana Road, Tambol Khanham,
 Amphur U-Thai
 Ayutthaya 13210, Thailand
 Tel: 66 35 33 4000
 Fax: 66 35 22 7088

Jakarta, Indonesia
 PT. Capsugel Indonesia
 Jl. Raya Bogor Km 42
 P.O. Box 15/CBI
 Cibinong 16916, Jawa Barat
 Indonesia
 Tel: 62-21-875-2226
 Fax: 62-21-875-2227

Suzhou, China
 Suzhou Capsugel Ltd.
 No. 369 Suhong Middle Road,
 Suzhou Industrial Park
 Jiangsu Province, P.R.China, 215021
 Tel: 86-512-6258-5188
 Fax: 86-512-6258-9188

Sales Office Hong Kong
 16/F, Stanhope House, 738 Kings Road
 North Point, Hong Kong
 Tel: (852) 2811 9711, 2562 5567 (direct)
 Fax: (852) 2214 1737

Sagamihara, Japan
 Capsugel Japan KK
 3-36 Minamihashimoto, 4-chome
 Sagamihara City
 Kanagawa Prefecture
 229-1133 Japan
 Tel: 81-427-00-6700
 Fax: 81-427-00-6719

Sydney, Australia
 Capsugel
 32-40 Cawarra Road
 Caringbah NSW, 2229
 Australia
 Tel: 61-2-9710-6642
 Fax: 61-2-9710-6283

For information on our other machines and services designed to speed your pharmaceutical development time, such as the Xcelodose™ Precision Micro-filling System, the CFS 1000 and the LEMS 60™, please visit us at www.capsugel.com.